COLD FUNCTION TEST GUN # A-7 SERIAL # 71126318

TRACKING # TLW 1216V DATE 9/10/03

- # 1 PLACE GUN A-7 AND 100 ROUND OF AMMO IN FREEZER FOR 15 HOURS AT -20 DEGREE F.
- # 2 TOOK GUN A-7 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO SHOOTING ROOM AND SHOT 20 ROUNDS LOADING ONE ROUND IN CHAMBER AND THREE ROUNDS IN THE MAGAZINE. CYCLED SAFETY EVERY FOURTH ROUNDS.

PUT GUN A-7 BACK IN FREEZER FOR 2 HOURS.

3 AFTER 2 HOURS TOOK GUN A-7 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM AND SHOT 20 ROUNDS.

LOADING ONE ROUND IN CHAMBER AND THREE ROUNDS IN THE MAGAZINE.

CYCLED SAFETY EVERY FOURTH ROUNDS.

PUT GUN A-7 BACK IN FREEZER FOR 2 HOURS.

4 AFTER 2 HOURS TOOK GUN A-7 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM AND SHOT 20 ROUNDS.

LOADING ONE ROUND IN CHAMBER AND THREE ROUNDS IN THE MAGAZINE.

CYCLED SAFETY EVERY FOURTH ROUNDS.

PUT GUN A-7 BACK IN FREEZER FOR 2 HOURS.

5 AFTER 2 HOURS TOOK GUN A-7 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOMS.

AND SHOT 20 ROUNDS.

LOADING ONE ROUND IN CHAMBER AND THREE ROUNDS IN THE MAGAZINE

CYCLED SAFETY EVERY FOURTH ROUNDS.

PUT GUN A-7 BACK IN FREEZER FOR 2 HOURS.

6 AFTER 2 HOURS TOOK GUN A-7 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM AND SHOT 20 ROUNDS.

LOADING ONE ROUND IN CHAMBER AND THREE ROUNDS IN THE MAGAZINE CYCLED SAFETY EVERY FOURTH ROUNDS

CLEANED AND INSPECTED (EVERYTHING OK)
MALFUNGTIONS