

CONFIDENTIAL

Remington Arms Company Inc.
RESEARCH & DEVELOPMENT TECHNICAL CENTER
315 WEST RING ROAD
ELIZABETHTOWN, KY 42701

TLW 1005

Table of Ammunition to use for Thermal Tests			
Manufacturer	Type	Code	RAC
Remington	High Velocity Game Load	HV12-6	20065
Remington	Premiere Target	ST812LH-8	20252
Remington	Gun Club Light Target	GC12L-8	20230
Remington	Express Mag. Buckshot	12SB00	20632
Remington	Heavy Field Load	RP12H-8	28120

Data Required:

- Shotgun serial number
- Cycle time for each test condition
- Temperature records throughout each cycle. Use the chart feature on the freezer and a temperature-recording device for the chamber.
- Malfunctions type and rates both pre- and post thermal cycles.
- Observations made on cracks, creep or other noteworthy items
- TLW Number
- Testers' Names

TLW1005V - Heat & Humidity Test:

This test evaluates the effects of high heat combined with high humidity such as might be found in some shooting environments. Of special interest for this test is the effect of heat and humidity on the wood stock.

Method:

- Store the gun and test ammunition for a minimum of six hours at a temperature of +100°F and 80-90% Relative Humidity.
- Test each firearm after removing from the chamber.
- Use ammunition schedule as listed in the table below.
- Fire 20 rounds of ammunition. Return the firearm to the chamber. Wait 2 hours and repeat procedure until all 100 rounds have been fired.
- Do not perform maintenance during the 100 round cycles.
- Cycle the safety from fire to safe every 5 rounds.
- After 100 rounds have been fired through the firearm, disassemble, thoroughly inspect, clean and lubricate.

Table of Ammunition to use for Thermal Tests			
Manufacturer	Type	Code	RAC
Remington	High Velocity Game Load	HV12-6	20065
Remington	Premiere Target	ST812LH-8	20252
Remington	Gun Club Light Target	GC12L-8	20230
Remington	Express Mag. Buckshot	12SB00	20632
Remington	Heavy Field Load	RP12H-8	28120

Data Required:

- Record temperature and exposure times
- Record all malfunctions
- Record damage noted during inspection
- Record Testers' Names
- Record TLW Number

J.R. Snedeker

Page 24 of 29

12/17/02

Remington Confidential

Revision #1.0

C:\Program Files\TCDI\CrackerLoaderREM\REME\email\rawblob\20060120133028A00024914.doc

Subject to Protective Order - Williams v. Remington

BARBER 5.30.06R00020845**ETE00001044**