

PROJECT # 241095
AMMO REMINGTON 30.06
180 GRAIN R30064

COLD FUNCTION TEST
GUN # C 22 SERIAL # 71001745

TRACKING # TLW 0300AH

AMMO REMINGTON 30.06 180 GRAIN R30064

LOT # K28GB3714

MALFUNCTIONS

MALFUNCTIONS

1 1
2 2
3 3
4 4
5 5 CYCLED SAFETY
6 6
7 7
8 8
9 9
10 10 CYCLED SAFETY
11 1
12 2
13 3
14 4
15 5 CYCLED SAFETY
16 6
17 7
18 8
19 9
20 10 FIRED 20 ROUNDS AT -20 DEGREE F.
21 1
22 2
23 3
24 4
25 5 CYCLED SAFETY
26 6
27 7
28 8
29 9
30 10 CYCLED SAFETY
31 1
32 2
33 3
34 4
35 5 CYCLED SAFETY
36 6
37 7
38 8
39 9
40 10 FIRED 20 ROUNDS AT -20 DEGREE F.
41 1
42 2
43 3
44 4
45 5 CYCLED SAFETY
46 6
47 7
48 8
49 9
50 10 CYCLED SAFETY

51 1
52 2
53 3
54 4
55 5 CYCLED SAFETY
56 6
57 7
58 8
59 9
60 10 FIRED 20 ROUNDS AT -20 DEGREE F.
61 1
62 2
63 3
64 4
65 5 CYCLED SAFETY
66 6
67 7
68 8
69 9
70 10 CYCLED SAFETY
71 1
72 2
73 3
74 4
75 5 CYCLED SAFETY
76 6
77 7
78 8
79 9
80 10 FIRED 20 ROUNDS AT -20 DEGREE F.
81 1
82 2
83 3
84 4
85 5 CYCLED SAFETY
86 6
87 7
88 8 MISFIRE
89 9
90 10 CYCLED SAFETY
91 1
92 2
93 3
94 4
95 5 CYCLED SAFETY
96 6
97 7
98 8
99 9
100 10 FIRED 20 ROUNDS AT -20 DEGREE F.

83

ET08363

PROJECT # 241095
AMMO REMINGTON 30.06
180 GRAIN R30064

COLD FUNCTION TEST
GUN # C 22 SERIAL # 71001745

TRACKING # TLW 0300AH

- # 1 PLACE GUN C 22 AND 100 ROUNDS OF AMMO IN FREEZER FOR 14 HOURS AT -20 DEGREE F.
- # 2 TOOK GUN C 22 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO SHOOTING ROOM AND SHOT 20 ROUNDS
LOADING ONE ROUND IN CHAMBER AND FOUR ROUNDS IN MAGAZINE.
CYCLED SAFETY EVERY FIVE ROUNDS.
PUT GUN C 22 BACK IN FREEZER FOR 2 HOURS.
- # 3 AFTER 2 HOURS TOOK GUN C 22 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM
AND SHOT 20 ROUNDS.
LOADING ONE ROUND IN CHAMBER AND FOUR ROUNDS IN MAGAZINE.
CYCLED SAFETY EVERY FIVE ROUNDS.
PUT GUN C 22 BACK IN FREEZER FOR 2 HOURS.
- # 4 AFTER 2 HOURS TOOK GUN C 22 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM
AND SHOT 20 ROUNDS.
LOADING ONE ROUND IN CHAMBER AND FOUR ROUNDS IN MAGAZINE.
CYCLED SAFETY EVERY FIVE ROUNDS.
PUT GUN C 22 BACK IN FREEZER FOR 2 HOURS.
- # 5 AFTER 2 HOURS TOOK GUN C 22 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM
AND SHOT 20 ROUNDS.
LOADING ONE ROUND IN CHAMBER AND FOUR ROUNDS IN MAGAZINE.
CYCLED SAFETY EVERY FIVE ROUNDS.
PUT GUN C 22 BACK IN FREEZER FOR 2 HOURS.
- # 6 AFTER 2 HOURS TOOK GUN C 22 AND 20 ROUNDS OF AMMO OUT OF FREEZER TO TO THE SHOOTING ROOM
AND SHOT 20 ROUNDS.
LOADING ONE ROUND IN CHAMBER AND FOUR ROUNDS IN MAGAZINE.
CYCLED SAFETY EVERY FIVE ROUNDS.
ROUND 88, SECOND ROUND OUT OF MAGAZINE (MISFIRE)
- # 7 CLEANED AND INSPECTED OK STEVE WADE

83

ET08364

GUN - C22
DATE - 12/12/00
SHOOTER - SW

M/710 100 ROUND

PROJECT# 241095

Cold Function Test

TLW03004H

AMMO TYPE - 30.06 R30064 180 Gr.		LOT# - K28G83714	
MALFUNCTIONS		MALFUNCTIONS	
1	1	51	3
2	2	52	4
3	3	53	1
4	4	54	2
5	1 ✓	55	3 ✓
6	2	56	4
7	3	57	1
8	4	58	2
9	1	59	3
10	2 ✓	60	4 ✓
11	3	61	1
12	4	62	2
13	1	63	3
14	2	64	4
15	3 ✓	65	1 ✓
16	4	66	2
17	1	67	3
18	2	68	4
19	3	69	1
20	4 ✓	70	2 ✓
21	1	71	3
22	2	72	4
23	3	73	1
24	4	74	2
25	1 ✓	75	3 ✓
26	2	76	4
27	3	77	1
28	4	78	2
29	1	79	3
30	2 ✓	80	4 ✓
31	3	81	1
32	4	82	2
33	1	83	3
34	2	84	4
35	3 ✓	85	1 ✓
36	4	86	2
37	1	87	3
38	2	88	4
39	3	89	1
40	4 ✓	90	2 ✓
41	1	91	3
42	2	92	4
43	3	93	1
44	4	94	2
45	1 ✓	95	3 ✓
46	2	96	4
47	3	97	1
48	4	98	2
49	1	99	3
50	2 ✓	100	4 ✓

Note: Magazine box hard to get out and put back in.